



Code	LT	X3	X1 min	X max	TJ	TH (x 2 pcs.)	ØC
CCLDL1DQT08029	290	143	>=38	<=214	SF151551000J	REQ50606501100U	5-6-7-8
CCLDL1DQT08034	340	168	>=38	<=264	SF151551000J	REQ50607501300U	5-6-7-8
CCLDL1DQT08039	390	193	>=38	<=314	SF151551000J	REQ50608501450U	5-6-7-8
CCLDL1DQT08044	440	218	>=38	<=364	SF151551000J	REQ50609501650U	5-6-7-8
CCLDL1DQT08049	490	243	>=38	<=414	SF151551000J	REQ50610501850U	5-6-7-8

Hot runner manifold with: No.8 aligned drops and double heater circuit (No.1 heater on injection bushing side, No.1 heater on nozzle side)

No.8 drops aligned - No. 1+1 heater
 TH = Armoured tubular heater 6 x 6 mm
 TJ = J-type thermocouple with M4 fixing device
 B = Min. 10 mm.(variable depending on the height of the nozzle head)
 NOTE : centres for nozzle head of Ø26