



Code	ØC	X3	X min	Y min	TJ	TH (x 2pcs.)	X * Y min ÷ max
CCLP1BQT06007	5 ÷ 10	TEC	54	27	SF151551000J	REQ506 SV W DIS+	÷ 7.225
CCLP1BQT06012	5 ÷ 10	TEC	54	27	SF151551000J	REQ506 SV W DIS+	7.226 ÷ 12.769
CCLP1BQT06022	5 ÷ 10	TEC	54	27	SF151551000J	REQ506 SV W DIS+	12.770 ÷ 21.904
CCLP1BQT06033	5 ÷ 10	TEC	54	27	SF151551000J	REQ506 SV W DIS+	21.905 ÷ 33.856

Hot runner manifold with: No.6 drops on two rows and double heater circuit (No.1 heater on injection bushing side, No.1 heater on nozzle side)

No.6 drops on two rows - No.1+1 heater
 TH = Armoured tubular heater 6 x 6 mm
 TJ = J-type thermocouple with M4 fixing device
 B = Min. 10 mm.(variable depending on the height of the nozzle head)
 NOTE : centres for nozzle head of Ø26
 TEC = Please ask EMP Technical Office